



4610 Littlejohn Street, Baldwin Park, CA 91706
(800) 321-5981, (626) 960-9822 Phone
(626) 960-8700 Fax
info@medinuclear.com
www.medinuclear.com

The Effect of Inertial Changes on Particle Size Distribution in Nebulizer Systems: The Impact on Clinical Aerosol Studies

*GA Francken, AD Waxman, R Potter and R King.
Cedars-Sinai Medical Center, L.A., CA, and Medi/Nuclear® Corp., Baldwin Park, CA.*

Previous studies have shown that aerosol particle size plays a major role in determining the quality of aerosol images produced from a nebulized radiopharmaceutical. The purpose of the current study was to compare a standard commercially available nebulizer with two modified systems with regard to particle size distribution and quality of aerosol image production.

Methods: A standard nebulizer system was modified to reduce the percent of particles greater than 1 micron by a series of directional changes which eliminated larger particles based upon inability of these particles to change direction rapidly. The unmodified nebulizer produced 35% of particles between 1.0-5.0 microns. The first modified system produced 2.93% of particles between 1.0-5.0 microns and the third produced 8.2% of particles between 1.0-5.0 microns.

43 consecutive patients were randomized to 1 of 3 nebulizers. Images were evaluated by 2 observers blinded as to the nebulizer used. A semi-quantitative rating system was employed for 3 characteristics of the images including overall image quality, peripheral penetration and central deposition.

In all three categories, the nebulizer demonstrating the lowest percentage of particles between 1.0-5.0 microns achieved the highest rating for all 3 categories. Central deposition was almost non-existent with this nebulizer.

We conclude that inertial modification of particle size to less than 1 micron is feasible. The reduction in particle size directly contributes to the improved quality of aerosol images resulting in increased peripheral penetration, decreased central deposition, and overall image quality.

Presented at the 22nd Annual Western Regional Meeting of the Society of Nuclear Medicine, Monterey, CA.

Dr. Francken was given the Young Investigators Award for the presentation.

September 18-21, 1997