

INSTA/VENT™ and UNIVERSAL MOBILE CABINETS

Medi/Nuclear's® proprietary Insta/Vent™ and Universal mobile cabinets are specifically designed for technologist safety and convenience. Each system is a complete point-of-care workstation that allows easy patient access and the convenient storage of both 'hot' and cold supplies and materials. A working countertop provides space for paperwork, doses, and/or patient belongings. Lead shield and oxygen tank holder sold separately.

INSTA/VENT MOBILE CABINET #IV-603

SPECIFICATIONS

Weight: 139 lbs.
Dimensions: 20 ¾ in. L x 15 ¼ in. W x 39 in. H
Decay Area: Full 1/8 in. lead lining

FEATURES

- ◆ Sliding shelf locks shield in place and brings it close to a patient.
- ◆ Large storage area for unused kits or supplies.
- ◆ Decay area with full lead lining.
- ◆ Secure drop down door.
- ◆ Veneer exterior for ease of cleaning.
- ◆ Swivel wheels.
- ◆ Locking front wheels.
- ◆ Flat top for managing paperwork and supplies.
- ◆ Oxygen tank holder (optional)



UNIVERSAL MOBILE CABINET™ #UV-204

SPECIFICATIONS

Weight: 110 lbs.
Dimensions: 19 1/2 in. L x 14 in. W x 26 in. H
Decay Area: Full 1/8 in. lead lining

FEATURES

- ◆ Storage area for unused kits and supplies.
- ◆ Locking decay area with full lead lining.
- ◆ Powder coating for ease of cleaning.
- ◆ Rubber bumpers.
- ◆ Swivel wheels.
- ◆ Locking front wheels.
- ◆ Flat top for managing materials.
- ◆ Swing Arm for Aero/Vent™ lead shields (optional)
- ◆ Oxygen tank holder (optional).



Model No.	ORDERING INFORMATION	Qty
IV-603	Insta/Vent™ Mobile Cabinet	Each
UV-204	Universal Mobile Cabinet™	Each
UV-204-SA	Universal Mobile Cabinet™ with Swing Arm	Each

www.medinuclear.com



Medi/Nuclear® Corporation
800.321.5981 / 626.960.9822 / 626.960.8700 Fax
info@medinuclear.com
4610 Littlejohn Street, Baldwin Park, CA 91706