



4610 Littlejohn Street, Baldwin Park, CA 91706  
(800) 321-5981, (626) 960-9822 Phone  
(626) 960-8700 Fax  
info@medinuclear.com  
www.medinuclear.com

Thank you for inquiring about the lifespan of Medi/Nuclear's® products.

As Medi/Nuclear's® radioaerosol delivery products, Xenon delivery products and aerosol drug delivery products do not require an expiration date, we are pleased to provide you with information regarding their shelf life.

The shelf life for Medi/Nuclear's® products is 7 years. However, as product components are made from high quality, long lasting materials, kits are generally safe to use for several years beyond that, when stored properly.

To enhance the lifespan of Medi/Nuclear's® radioaerosol kits, store them in moderate room temperatures, avoiding extreme heat, cold or light. As with any product, best practices recommend rotating stock so oldest products are used first.

To learn more about Medi/Nuclear's® products you are invited to visit the Global Unique Device Identification Database (GUDID) website at <https://accessgudid.nlm.nih.gov/>. There you will find information about device characteristics, materials, storage and more. Should additional product information be needed, please let us know. As always, we're here to be of assistance.

Thank you once again,

A handwritten signature in blue ink, appearing to read "Suzanne Hammersmith".

Suzanne Hammersmith, CEO  
Medi/Nuclear® Corporation